

# Squid<sup>TM</sup>



## The New Liquid Embolic Device

Squid is a non adhesive liquid embolic agent for embolization of brain Arteriovenous Malformations (AVM). It is composed of EVOH (ethylene vinyl alcohol copolymer) with suspended micronized Tantalum powder for radiopacity and DMSO (dimethyl sulfoxide) solvent. Squid must be injected through a DMSO compatible microcatheter.


**SQUID EMBOLIZATION RANGE :  
THE POWER OF VISIBILITY, THE CHOICE OF FLUIDNESS**

**emboflu<sup>AG</sup>**

Route des Avouillons 30 ♦ 1196 GLAND ♦ SWITZERLAND  
Tel: +41 22 354 07 10 ♦ info@emboflu.ch ♦ www.emboflu.ch

## A UNIQUE MICRONIZATION PROCESS

EMBOFLU developed a specific micronization process to minimize Tantalum powder grain size in the SQUID solution. This enables a SLOWER precipitation of the radiopaque powder which stays LONGER in the SQUID suspension.



## BENEFITS of the SQUID unique micronization process

- ◆ High HOMOGENEITY of SQUID suspension - eliminates aggregates formation which can cause microcatheter blockage.
- ◆ HOMOGENIC radiopacity - reduces discrepancy between saturated radiopaque zones and non-saturated zones, for improved assessment of AVM angioarchitecture.
- ◆ High STABILITY over time - allows improved VISIBILITY for longer INJECTION TIMES.


# 4 SQUID FORMULAS

## STANDARD VISCOSITY


**SQUID18:** Standard version of SQUID  
For standard AVM embolization.

**SQUID18LD:** Low Density  
30% less radiopacity than standard SQUID18.

- ◆ For a better assessment of the AVM vasculature and the volume of embolic liquid injected.
- ◆ To avoid the over saturated radiopaque injected zones "Flash effect".



Courtesy of Dr Gal  
Odense, Denmark



SQUID 12  
reduced  
viscosity enables  
penetration into  
microvessels and  
nidus of AVM

SONIC  
microcatheter

Courtesy of Dr. Amsalem  
Haifa Israel

## LOW VISCOSITY

**SQUID12:** With a lower viscosity,  
this version is more fluid than the standard  
formula to allow:

- ◆ Deeper penetration into the nidus.
- ◆ To reach distal microvessels and inject through small feeders.

**SQUID12LD:** Low Density

- ◆ This formula has the same viscosity as SQUID 12, but 30% less Tantalum.

- ◆ To enable assessment of the AVM vasculature and the volume of liquid embolic injected.



Squid 12LD enables the visibility  
of the sonic distal marker

# SQUID<sup>TM</sup>

## The New Liquid Embolic Device

### ORDERING INFORMATION:

REFERENCES	DESCRIPTION	CONTENT OF EACH REFERENCE
SQUID 18	STANDARD VISCOSITY SQUID	<ul style="list-style-type: none"><li>◆ One 1.5ml vial of SQUID</li><li>◆ One 1.5 ml vial of DMSO</li><li>◆ One 1cc Blue syringe for DMSO</li><li>◆ Two 1cc White syringes for SQUID</li><li>◆ Two syringe adapters</li></ul>
SQUID18LD	STANDARD VISCOSITY and LOW DENSITY SQUID: 30% less radiopacity than SQUID 18	
SQUID 12	LOW VISCOSITY SQUID	
SQUID12LD	LOW VISCOSITY and LOW DENSITY SQUID: 30% less radiopacity than SQUID 12	

**SQUID EMBOLIZATION RANGE :  
THE POWER OF VISIBILITY, THE CHOICE OF FLUIDNESS**



**emboflu<sup>AG</sup>**

Route des Avouillons 30 ◆ 1196 GLAND ◆ SWITZERLAND  
Tel: +41 22 354 07 10 ◆ info@emboflu.ch ◆ www.emboflu.ch