

goldballoon

balt's exclusive detachable balloon




avm* treatment

goldballoon

Detachable balloon

Detachable latex balloons designed for the embolization of arterio-venous malformations in the neurological & peripheral regions



with the courtesy of Dr Rujimas Khumtong

Unique design

Latex compliant balloons with a RX marker ensuring an accurate placement & conformability

Latex valve enabling quick fitting on catheters & highly resistant to avoid leakage & deflation

Unique delivery setting

The only detachable balloon on the market compatible with the following microcatheters:

MABDTE
BALTACCIBDPE
BALTACCIBDPE 100

ordering information

A variety of diameters available

Reference	Volume (mL)	Inflated balloon ØxL (mm)	Collar Ø (mm)	Minimum guiding catheter I.D.
GOLDBAL1	0,25	7 x 10	0,90	1,75mm (0.07")
GOLDBAL2	0,60	7 x 20	0,90	
GOLDBAL3	0,50	9 x 11	0,90	2,20mm (0.09")
GOLDBAL4	0,75	9 x 16	0,90	
GOLDBAL5	2,50	12 x 28	0,90	2,80mm (0.11")
GOLDBAL6	3,00	15 x 25	1,05	3,10mm (0.12")

Delivery catheters

Reference	O.D. (F) prox / dist	Total length (cm)
MABDTE	2,7 / 1,8	165
BALTACCIBDPE	2,7 / 2,4	163
BALTACCIBDPE100	2,7 / 2,4	100

* Arteriovenous Malformation

GOLDBALLOON is designed for the embolisation of arterio-venous malformations in the neurological and peripheral regions. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/EEC amended by 2007/47/EC). Manufactured by BALT EXTRUSION SAS. Carefully read the instructions for use before use. Not reimbursed. First CE marking:2003. MABDTE and BALTACCIBDPE are catheters that allow the introduction of detachable balloons (GOLDBALLOON) and are indicated for use in the neurovasculature and peripheral system. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/EEC amended by 2007/47/EC). Manufactured by BALT EXTRUSION SAS. Carefully read the instructions for use before use. First CE marking:2007. The content of this document, in particular data, information, trademarks and logos is BALT SAS and affiliates' sole property. © 2020 BALT SAS and affiliates, all rights reserved. All representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT SAS and affiliates' copyrights and other intellectual proprietary rights. This document with associated pictures is non-contractual and is solely dedicated to healthcare professionals and BALT's distributors (BALT's supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. DC048GB (06/20).

Balt

10, rue de la Croix Vigneron, 95160 Montmorency, France

Tél. : +33 (0)1 39 89 46 41

Fax : +33 (0)1 34 17 03 46

www.balt-corp.com

