

ANEURYSMS CAN BE COMPLEX

SIZING
SHOULD BE
SIMPLE

CONTOUR

Neurovascular
System™



Cerus
Endovascular



CONTOUR
 Neurovascular
 System™
5
 SKUs
 2.0 - 10.5mm

V

Competitive
 Intrasaccular
 Device
44
 SKUs
 3.0 - 11.0mm

The simple sizing of the Contour Neurovascular System™ requires only the measurement of the maximum aneurysm neck and diameter – this means, you are able to select from only 5 Contour devices which cover the same full size range as 44 alternative competitive intra-saccular devices.

Contour Neurovascular System™ Sizing Chart

Product Code	MC	Device Diameter (mm)	Aneurysm Neck (mm)	Aneurysm Width (mm)
CNS005-15	0.027"	5.0	2.0 - 3.0	2.0 - 3.5
CNS007-15	0.027"	7.0	3.0 - 5.0	3.0 - 5.5
CNS009-15	0.027"	9.0	4.0 - 6.0	5.0 - 7.5
CNS011-15	0.027"	11.0	5.0 - 8.0	7.0 - 8.5
CNS014-15*	0.027"	14.0	7.0 - 10.0	8.0 - 10.5

*Coming soon

Contour Neurovascular System™ is a registered trademark of Cerus Endovascular Ltd. Prior to use, refer to the Instructions for Use supplied with products for indications, contraindications, warnings, precautions, and suggested procedure.

The Contour Neurovascular System™ meets the essential requirements referred to in Article 3 of the MDD 93/42/EEC and are labelled with the CE Marking of conformity. Product not for sale in the United States. ©2020 Cerus Endovascular, Inc. All rights reserved. APM014-01.A

Cerusendo.com



Cerus
 Endovascular