

WHEN EVERY SECOND COUNTS


NEQSTENT™

Coil Assisted
Flow Diverter


SIZING
SHOULD BE
SIMPLE



Cerus
Endovascular



- NQS414-XX*
- NQS21411-15
- NQS21409-15
- NQS21407-15



When dealing with complex bifurcation aneurysms, **every second counts**. Neqstent provides quick and simple sizing, by only needing to measure the neck of the aneurysm.

With just 4 sizes, you can treat aneurysm necks from 3.0 to 10.0 mm.

Neqstent Coil Assisted Flow Diverter Sizing Chart

Product Code	MC	Device Diameter (mm)	Aneurysm Neck (mm)
NQS21407-15	0.021"	7.0	3.0 - 5.0
NQS21409-15	0.021"	9.0	4.0 - 6.0
NQS21411-15	0.021"	11.0	5.0 - 8.0
NQS414-XX*	0.027"	14.0	7.0 - 10.0

**Coming soon*

Prior to use, refer to the Instructions for Use supplied with products for indications, contraindications, warnings, precautions, and suggested procedure. Neqstent meets the essential requirements referred to in

Article 3 of the MDD 93/42/EEC and are labelled with the CE Marking of conformity. Product not available for sale in the United States. ©2020 Cerus Endovascular, Inc. All rights reserved. APM403-01-A

Cerusendo.com



Cerus Endovascular