



Leading BGC Compatibility

Market Leading Intermediate & Aspiration Catheters	OD in.	OD mm	ID in.	ID mm	Q'Apel Walrus 0.087" ID	Stryker Flowgate 2.084" ID
Penumbra Red 62	0.076	1.93	0.062	1.57	✓	✓
Penumbra Red 68	0.084	2.13	0.068	1.73	✓	✗
Penumbra Red 72	0.085	2.16	0.072	1.83	✓	✗
Penumbra Jet 7	0.085	2.16	0.072	1.83	✓	✗
Penumbra Jet D Reperfusion Catheter	0.065	1.65	0.064 P/ 0.054 D	1.63 P/ 1.37 D	✓	✗
Penumbra ACE 060	0.071 P/ 0.068 D	1.80 P/ 1.73 D	0.060	1.52	✓	✓
Penumbra ACE 064	0.084 P/ 0.080 D	2.13 P/ 2.03 D	0.068 P/ 0.064 D	1.73 P/ 1.63 D	✓	✓
Penumbra ACE 068	0.084	2.13	0.068	1.73	✓	✓
Stryker AXS Catalyst 5	0.073 P/ 0.069 D	1.85 P/ 1.73 D	0.058	1.47	✓	✓
Stryker AXS Catalyst 6	0.079 P/ 0.071 D	2.01 P/ 1.80 D	0.060	1.52	✓	✓
Stryker AXS Vecta	0.085 P/ 0.082 D	2.16 P/ 2.08 D	0.071	1.8	✓	✗
Stryker Vecta 74	0.087	2.21	0.074	1.88	✗	✗
Medtronic React 68	0.083	2.11	0.068	1.73	✓	✓
Medtronic React 71	0.0855	2.17	0.071	1.8	✓	✗
Microvention Sofia 5F	0.068	1.73	0.055	1.4	✓	✓
Microvention Sofia Plus 6F	0.083 P/ 0.082 D	2.11 P/ 2.08 D	0.070	1.78	✓	✗
Microvention Sofia Flow Plus	0.0825 P/ 0.0815 D	2.10 P/ 2.07 D	0.070	1.78	✓	✓
Imperative Care Zoom 71	0.083	2.1	0.071	1.8	✓	✓
Imperative Care Zoom 88	0.110 P/ 0.106 D	2.79 P/ 0.106 D	0.088	2.24	✗	✗
Cerenovous Large Bore Catheter LBC	0.0825 P/ 0.081 D	2.10 P/ 2.06 D	0.071	1.8	✓	✓

Walrus is compatible with all carotid stent systems.

REF

BG8087-090

BG8087-095

We're Q'Apel, a medical technology company that creates solutions. More specifically, we design novel access device technology for vascular interventions and unmet clinical needs. Because in the precious seconds that surround a stroke emergency, clinicians need technology that delivers. That's where we come in.

Get in touch | info@qapelmedical.com | [linkedin.com/company/q-apel-medical](https://www.linkedin.com/company/q-apel-medical)

Q'Apel Medical

4245 Technology Drive

Fremont, CA 94538

510.738.6255

[qapelmedical.com](https://www.qapelmedical.com)