

KRONOS™

ELECTROCAUTERY DEVICE



<https://app.nexposai.com/item/kronos/>

- ▶ Easy to Use
- ▶ Rapid Bleeding Control
- ▶ Portable and Disposable



The only single-use electrocautery device that
can cauterize deep tissue through a guide needle


singlepass®

Kronos Electrocautery Device

Ordering Information

Part Number	Model	Working Length	Outer Diameter (OD)	Description
900001SP-01	K10	17cm	0.040in	Kronos Electrocautery Device; Compatible with 10cm guide needle
900001SP-02	K15	22cm	0.040in	Kronos Electrocautery Device; Compatible with 15cm guide needle

Device compatible with up to 17-gauge biopsy guide needle up to 15cm long.



The Single Pass Kronos Electrocautery device is indicated for solid organ (i.e., liver, kidney, lungs, etc.) biopsy procedures. The Single Pass Kronos Electrocautery device is intended to be used in biopsy and minor invasive procedures to control bleeding by use of heated probe tip that is applied directly to the target tissue.



Phone: 1-888-311-6777 Email: info@singlepass.co Web: singlepass.co Address: Single Pass, Inc. 105 N Pointe Dr, Suite B Lake Forest, CA 92630 USA

MKTG-0001SP Rev C550903MDR201