

VIRTUAL PORTS

Giving you room to maneuver



Introducing the First End to End Solution for Organ Maneuvering during Laparoscopic Surgery

Traditional Laparoscopy | Robotic Assisted Laparoscopy | Single Port Laparoscopy

EndoLift


**Ultra-versatile organ and tissue maneuvering system
based on proprietary Micro Anchoring™ Technology**

The **EndoLift** system is a single hand, organ suspension, lifting, and retracting device for Traditional, Robotic Assisted or Single Port Laparoscopic surgery. Applied through a 5mm port, **EndoLift** is latched on to the peritoneum retro to the organ designated for repositioning, allowing organ retracting, lifting or suspension at the desired angle and location. The **EndoLift** replaces multiple devices used in laparoscopy procedures today and reduces surgeon dependency on assistants, as well as streamlining procedures and freeing the laparoscopy access port for other operational needs.

EndoLift's unique design enables a flexible choice of the direction and angle of lifting, doesn't require anchoring to any external device, which may limit the patient's position and adds unwanted operating table clutter.



Two anchoring clips to
attach to the internal cavity
wall (peritoneum)




Interconnected by a
telescopic tube to support
organs with no trauma




Single hand
positioning

Multi-port Laparoscopy

Virtual Ports Freedom To Operate



High Surgeon Dependency on
Assistance and Tools




One Surgeon | One Hand | One Tool

Multiple Clinical Applications


Advanced capabilities to maneuver intraoperatively aiming to reduce operating clutter, cleared use in a variety of clinical applications, providing greater flexibility and agility for a more accurate, safer outcome.



Upper GI




Colon & Rectal



HPB



GYN



Urology


THE FREEDOM TO OPERATE REVOLUTION



All in one multiple
application device



Unmatched
organ and tissue
maneuvering



Keeping ports free
and available



Independent
single hand
application



Reduced cost and
complications

MAKING AN IMPACT

Our solutions support better patient care through shorter, more successful surgical procedures, less hospitalizations and adverse events and higher operations capacity a day.

avg

15%

Reduced OR time

up to

20%



Cost reduction

>30

Applications

>25K

Safely performed
surgeries

Ordering Information	Product Description	Product Code	Quantity
	EndoGrab Maneuvering Device - Sterile (pack of 10)	EG000A10	10
	EndoLift Suspension Device - Sterile (pack of 5)	EL000A5	5
	Applier 30cm for VP Lap Devices	AP030A	1
	Applier 45cm for VP Lap Devices	AP045A	1
	Applier PL 30cm for VP Lap Devices	APL030A	1
	Applier PL 45cm for VP Lap Devices	APL045A	1
	Protective Tube for Applier (pack of 2)	AP139A2	2

We believe that our End to End Solution for Organ Maneuvering during Laparoscopic Surgery adds value for patients and healthcare systems. Thank you for reviewing this information.

MKT-EL-01Rev.1EndoLift Brochure